

RELEASE OF DEVICE

by the then Governor of Jharkhand, Shri CP Radhakrishnan

PATENT

Indian Patent Office Application No. 202331062688.

ECONOMICS

Bill of materials : Rs. 15000 (Approx.) including parts and labour

Market Price : Rs. 25000 (Approx.)

Profit per unit : Rs. 10000 (Approx.)

Potential users: Lac processing centres, Quality Control & Assurance Labs, Food Regulatory Authorities, the Lac-based industry, the Paint Industry, and ICAR- KVK.

Technology Transferred: M/s Mansa Enterprises, Chandigarh

Authors

Nandkishore Thombare, Saurav Kumar, AP Bhondekar, Anupma Sharma, Ritesh Kumar, Sudeshna Bagchi, MF Ansari & Abhijit Kar

Published by

Director, ICAR- NISA, Namkum, Ranchi

(Sept., 2025)

Developed by ICAR NISA, Ranchi and CSIR CSIO, Chandigarh, under NABARD sponsored project

Contact For further Information

Director

ICAR- National Institute of Secondary Agriculture
Namkum, Ranchi 834 010

Email- director.nisa@icar.gov.in, director.icar.nisa@gmail.com Website: https://nisa.res.in/

LAC COLOR INDEX ANALYSER

An electronic, handheld, low cost and Make in India device for the quality analysis of lac resin





ICAR- NATIONAL INSTITUTE OF SECONDARY AGRICULTURE
Namkum, Ranchi- 834 010, Jharkhand, India




97. (0825) Anapuma Press and Process 집을 628799<mark>9194, e-mail : appranchi.1969@gma</mark>

LAC COLOR INDEX ANALYSER (LACILYSER)

WHAT IS COLOR INDEX?



The color index of lac is a standardized numerical or categorical value that represents the color characteristics of lac. It provides a simplified way to assess and communicate the color properties of lac, allowing for easy comparison and specification compliance in various industries.



HOW IS IT MEASURED CONVENTIONALLY?

According to the Indian Standard IS 6921:1973 (Methods of sampling and test for lac and lac products), the colour index of lac is estimated by visually comparing the lac solution with a standard iodine solution of known colour intensity. The lac solution appears darker, and it is diluted step by step with alcohol until its colour matches the iodine standard. The final colour index value is calculated based on the volume of alcohol added to achieve the match. The darker the sample, the higher its colour index will be.

LIMITATIONS OF CONVENTIONAL METHODS

- Time-consuming
- Subjective in analysis
- Required skilled workforce
- Accuracy and precision are compromised



HOW IS THE LACILYSER UNIQUE?

LaCilyser: The Lac Color Index Analyser is an integrated device designed to measure the color index of lac precisely. With its pre-calibrated system and resolution of color index ±1.0, the device offers accurate results across a wide range of lac color index values. Operating on a user-friendly interface, the device incorporates data acquisition and processing capabilities, providing quick and reliable measurements. The standardized preparation of the lac solution, as per IS:6921-1973, ensures consistency and adherence to recognized protocols. The LaCilyser is a valuable tool for industries requiring precise color values, contributing to quality control.

HOW ARE SAMPLES TESTED USING LACILYSER?

Step 1: Prepare lac samples as per BIS standards and pour in the sampling tube.

Step 2: Insert the sampling test tube in the testing chamber.

Step 3: Press the measuring button on the left side of the device. Color Index value is displayed on the screen.

SALIENT FEATURES

· No Titration Rapid Analysis

Minimal Safety Requirement

Precise and Accurate

Low Cost of Analysis

Low Maintenance

Low Operations Cost

Easy Sample Preparation

Easy Operation

Portable and Ergonomic

Customizable to any resins

SPECIFICATIONS

Measurement: Range 0 to 40 color index values

Analysis time : 5 sec

Accuracy : +/- 1.0 Color index
Precision : +/- 0.5 Color index
Resolution : - 0.1 Color index
Error : 1.0 Color index
Sensor Type : Spectral Sensor

Display : 2.5" LCD

Flexibility : Logically Programmable,

Customized and User friendly

Sensitivity : +/- 1.0 Color index

Repeatability: 0.01 decimal place value

Battery backup: 10 hours per full charge